

Contec® Healthcare TB1-3300™ Disinfectant

Ethyl Alcohol/Quat Cleaner Disinfectant

Contec® Healthcare TB1-3300™ Disinfectant is a blend of quaternary ammonium and ethyl alcohol available in 32 oz. (0.9L) and 1 gallon (3.8L) containers. TB1-3300 Disinfectant is available sterile and nonsterile. It is a ready to use cleaner/disinfectant that generates low residue.

TB1-3300 Disinfectant has a 1-minute bactericidal, virucidal, tuberculocidal and fungicidal contact time. It is bleach-free, making it an excellent choice for operators in critical environments, sterile compounding facilities, retail pharmacies, and laboratories.

A push-pull cap is included for controlled application. Contec Healthcare TB1-3300 Disinfectant is ideal for cleaning and disinfecting environmental surfaces in pharmacy cleanrooms and laboratory areas including work stations, isolators, Primary Engineering Controls (PECs), Secondary Engineering Controls (SECs), and other hard non-porous environmental surfaces.

Please contact your sales representative for further information.



| Features | Benefits |
|-------------------------------|--|
| One-step cleaner/disinfectant | <ul style="list-style-type: none"> • Compatible with common cleanroom surfaces • Great disinfectant/germicidal detergent for daily use |
| 1-minute contact time | <ul style="list-style-type: none"> • Easy to maintain compliance and appropriate dwell time |
| Low residue | <ul style="list-style-type: none"> • Less residue left behind after use • Residues are easily removed with common solutions* |
| 0.2 micron filtered | <ul style="list-style-type: none"> • Appropriate for use in critical environments • No filtering required before use |
| Push-pull cap | <ul style="list-style-type: none"> • Better application control of disinfectant |
| Linear tear packaging | <ul style="list-style-type: none"> • Easy to open even when wearing gloves • Facilitates transfer into more critical areas |
| Bleach-free | <ul style="list-style-type: none"> • Compatible with common controlled environment surfaces |

*Recommended rinse agents are isopropanol, ethanol, or sterile water.

| Product Specifications | |
|--|--|
| Odor | Alcohol |
| EPA-Registered | Yes |
| % Ethyl Alcohol | 72.5% |
| % Didecyl Dimethyl Ammonium Chloride | 0.33% |
| Detergency | Excellent |
| Density | 0.86* |
| pH | 3.7 at 20°C |
| Shelf Life Unopened (sterile and nonsterile) Opened (sterile and nonsterile) | 2 years 6 months (chemical potency) |
| HMIS Hazard Rating: | |
| Health | 2 |
| Flammability | 3 |
| Reactivity | 0 |
| Personal Protection | X |

*See SDS for more detailed information.

| Part No. | Description | Packaging |
|------------|--|-----------------|
| HCTB3300 | Contec Healthcare TB1-3300 Disinfectant, 32 oz. (0.9L), single bagged * | 12 bottles/case |
| HCTB3300IR | Contec Healthcare TB1-3300 Disinfectant, 32 oz. (0.9L), sterile, double bagged * | 12 bottles/case |
| HCTB3305 | Contec Healthcare TB1-3300 Disinfectant, 1 gallon (3.8L), single bagged | 4 bottles/case |
| HCTB3305IR | Contec Healthcare TB1 - 3300 Disinfectant, 1 gallon (3.8L), sterile, double bagged | 4 bottles/case |

Note: All units contain a push-pull cap.



Contec® Healthcare TB1-3300™ Disinfectant Organism Efficacy

Disinfectant Claims Summary

Bactericidal Claim – 1 minute

| |
|------------------------------------|
| Bordetella pertussis |
| Burkholderia cepacia |
| Campylobacter jejuni |
| Enterobacter aerogenes |
| Escherichia coli O157:H7 (E. coli) |
| Klebsiella pneumoniae |
| Legionella pneumophila |
| Listeria monocytogenes |
| Proteus vulgaris |
| Pseudomonas aeruginosa |
| Salmonella enterica |
| Shigella dysenteriae |
| Staphylococcus aureus |
| Streptococcus pyogenes (Strep) |
| Streptococcus pneumoniae |

Drug Resistant and Multi-Drug Resistant (MDR) Bacteria

| |
|---|
| Acinetobacter baumannii – Multi-Drug resistant (MDR) |
| Carbapenem resistant Enterobacter cloacae (CRE) (KPC) (MDR) |
| Carbapenem resistant Enterobacter cloacae (CRE) (New Delhi) (MDR) |
| Enterococcus faecium – Vancomycin resistant (VRE) |
| Escherichia coli – ESBL Producing (MDR) |
| Klebsiella pneumoniae – ESBL Producing (MDR) |
| Staphylococcus aureus – Methicillin resistant (MRSA) |
| Staphylococcus aureus – Vancomycin resistant (VRSA) |
| Staphylococcus aureus Panel – Vancomycin intermediate (VISA) |

Fungicidal Claim – 1 minute

| |
|----------------------------|
| Candida albicans |
| Trichophyton interdigitale |

Virucidal Claim – 1 minute

| |
|--------------------------------------|
| Adenovirus type 5 |
| Avian Influenza A Virus |
| Cytomegalovirus |
| Enterovirus EV-D68 |
| Hepatitis B Virus (HBV) |
| Hepatitis C Virus (HCV) |
| Herpes Simplex Virus Type 1 |
| Human Coronavirus |
| Human Immunodeficiency Virus (HIV-1) |
| Influenza Virus type A2 |
| Influenza B |
| Norovirus |
| Poliovirus serotype 1 (PV-1) |
| Respiratory Syncytial Virus |
| Rhinovirus |
| Rotavirus Strain WA |
| Vaccinia Virus |

Animal Pathogens

| |
|------------------------------|
| Avian Influenza A Virus |
| Bordetella bronchiseptica |
| Bovine Diarrhea Virus (BVDV) |
| Canine distemper Virus |
| Canine parvovirus |
| Duck Hepatitis B Virus |
| Feline calicivirus |
| Feline leukemia Virus |

Tuberculocidal Claim – 1 minute

| |
|---------------------|
| Mycobacterium bovis |
|---------------------|