

PROSAT® Sterile™ Low Endotoxin Wipes

Polynit Heatseal wipes with 70% USP grade IPA, 30% Water for Injection

PROSAT Sterile Low Endotoxin Wipes are made from 100% knitted polyester Polynit Heatseal wipes presaturated with a blend of 70% IPA and Water for Injection (WFI). Each batch of wipes has a guaranteed endotoxin level of less than <0.25EU/mL per device. Suitable for use in aseptic suites.

The sealed edge wipe produces a very low level of particles and fibers and is ideal for use near product contact surfaces.

Each pouch contains a small number of wipes eliminating any waste at the end of a session. Triple bagged for easy entry into the sterile suites. Validated sterile, these Low Endotoxin wipes are suitable for use in the most critical pharmaceutical cleanrooms.



FEATURE	BENEFIT
Cleanroom laundered 100% knitted interlock polyester wipes	<ul style="list-style-type: none"> • Very low levels of particles and fibers • Good abrasion and chemical resistance
Laser sealed edges	<ul style="list-style-type: none"> • Reduces fiber generation from the edge of the wipe
Presaturated with 70% IPA/WFI	<ul style="list-style-type: none"> • Controlled saturation levels greatly reduce solvent usage, VOC's and waste
Guaranteed endotoxin levels of <0.25EU/mL	<ul style="list-style-type: none"> • Suitable for use near product contact surfaces
Validated sterile via gamma irradiation to a 10 ⁻⁶ Sterility Assurance Level	<ul style="list-style-type: none"> • Suitable for use in ISO Class 3-4 cleanrooms
Small number of wipes per resealable pouch	<ul style="list-style-type: none"> • Low number of wipes eliminates waste • Pouch can be easily resealed throughout use, maintaining saturation levels



Part #	Description	Packaging
PSWE0002	PROSAT Sterile with 70% IPA/30% WFI, Polynit Heatseal 12" x12" (30.5 x 30.5cm)	30 wipes per pouch



P.O. Box 530 • Spartanburg, SC 29304 USA
Tel: +1-864-503-8333
Toll Free (US, Canada, PR): 1-800-289-5762
www.contecinc.com • wipers@contecinc.com