

PROSAT® Wipes

Presaturated nonwoven wipes

Contec® Healthcare PROSAT® wipes are low-lint wipes containing 70% USP grade isopropyl alcohol and 30% DI water. These nonwoven wipes are suitable for ISO 7 and 8 environments. Manufactured in Contec's ISO registered facility, these wipes are designed to remove the guesswork and inconsistency associated with spray and wipe processes making it possible to create a cleaning protocol that is repeatable, consistent, effective, and easier to implement, document, and validate.

Contec's PROSAT wipes are packaged in individually flow-wrapped, resealable pouches, and tamper-evident labels ensure product integrity and cleanliness, reducing solvent storage issues. These wipes are ideal for wiping supplies prior to pass through into the cleanroom and are convenient for wiping sterile gloves. In addition, PROSAT wipes provide residue removal on work surfaces outside the primary engineering control, after applying an EPA-registered one-step disinfectant cleaner. PROSAT wipes also allow for convenient cleanup following surface sampling.



Packaging facilitates easy opening and resealing to avoid dried wipes and potential contamination.

To learn more about Contec Healthcare's PROSAT Wipes for sterile compounding pharmacies, contact your Contec Healthcare sales representative or visit our website: www.contechealthcare.com.

Features	Benefits
Presaturated Wipes	<ul style="list-style-type: none"> • Contains consistent amount of 70% USP grade isopropyl alcohol (IPA) • Promotes process repeatability • Reduces VOC emissions often associated with spray options
Low-lint Construction	<ul style="list-style-type: none"> • Low-lint, meltblown polypropylene (MBPP) allows for even distribution of solution to intended surfaces • Low-lint, polyester/ cellulose blend fabric offers a stronger fabric, ideal for aggressive cleaning and wiping of rougher surfaces
Optimal Packaging	<ul style="list-style-type: none"> • When properly sealed, peel and reseal packaging eliminates evaporation concern during storage • Individually flow-wrapped, resealable pouches with tamper-evident labels reduce potential risk of contamination introduction into critical environments • Reduces solvent storage needs and potential spill hazards

Part No.	Description	Size	Packaging
HCPS2328	Presaturated meltblown polypropylene wipes with 70% IPA/30% DI water	9x11" (23x28cm)	30 wipes/pouch, 1 pouch/inner bag, 9 inner bags/outer bag, 3 outer bags/case; 810 wipes/case
HCPS0002	Presaturated cellulose/polyester wipes with 70% IPA/30% DI water	9x11" (23x28cm)	30 wipes/pouch, 1 pouch/inner bag, 9 inner bags/outer bag, 3 outer bags/case; 810 wipes/case