

# PROSAT® Sterile™ Low Endotoxin Wipes

Polynit Heatseal wipes with 70% USP grade IPA,  
30% Water for Injection

PROSAT Sterile Low Endotoxin Wipes are made 100% knitted polyester Polynit Heatseal wipes presaturated with a blend of 70% IPA and Water for Injection (WFI). Each batch of wipes has a guaranteed endotoxin level of less than <1EU/device. Suitable for use in aseptic suites.

The sealed edge wipe produces a very low level of particles and fibers and is ideal for use near product contact surfaces. Hanging slot in pouch makes it easy to use in isolators.

Each pouch contains a small number of wipes eliminating any waste at the end of a session. Double or triple bagged for easy entry into the sterile suites. Validated sterile, these Low Endotoxin wipes are suitable for use in the most critical pharmaceutical cleanrooms.



Features	Benefits
<b>Cleanroom laundered 100% knitted polyester wipes</b>	<ul style="list-style-type: none"> <li>• Very low levels of particles and fibers</li> <li>• Good abrasion and chemical resistance</li> </ul>
<b>Laser sealed edges</b>	<ul style="list-style-type: none"> <li>• Reduces fiber generation from the edge of the wipe</li> </ul>
<b>Presaturated with 70% IPA/WFI</b>	<ul style="list-style-type: none"> <li>• Controlled saturation levels greatly reduce solvent usage, VOC's and waste</li> </ul>
<b>Guaranteed endotoxin levels of &lt;1EU/device</b>	<ul style="list-style-type: none"> <li>• Suitable for use near product contact surfaces</li> </ul>
<b>Validated sterile via gamma irradiation to a 10<sup>-6</sup> Sterility Assurance Level</b>	<ul style="list-style-type: none"> <li>• Suitable for use in Grade A and B cleanrooms</li> </ul>
<b>Small number of wipes per resealable pouch</b>	<ul style="list-style-type: none"> <li>• Low number of wipes eliminates waste</li> <li>• Pouch can be easily resealed throughout use, maintaining saturation level</li> </ul>
<b>Double or triple bagged in linear tear packaging</b>	<ul style="list-style-type: none"> <li>• Each bag is easy to open even when wearing gloves</li> <li>• Facilitates cleaning products entering isolators and RABS</li> </ul>

Part No.	Description	Quantity
PSWE0001	PROSAT Sterile Polynit Heatseal with 70% IPA/30%WFI, 9" x 9" (23 x 23cm), triple bagged	10 wipes/pouch; 55 pouches/case