

# Contec® Critical Site™ Wipes

Sterile presaturated wipes designed for convenience in the pharmacy

Contec® Critical Site™ Wipes are 100% polyester 4" x 4" (10.2 x 10.2cm) wipes that contain 70% USP Grade Isopropanol and 30% USP Grade purified water. The low linting, nonwoven wipes are suitable for ISO Class 5 environments, are USP <797> compliant, and are exceptionally clean making them ideal for the most critical applications. Manufactured in Contec's ISO registered facility, the wipes are Validated Sterile per the Association for the Advancement of Medical Instrumentation (AAMI) 11137 Guidelines to a 10<sup>-6</sup> Sterility Assurance Level (SAL). Sterility documentation is available upon request.

Presaturated Contec Critical Site Wipes create a repeatable, consistent technique that allows pharmacists and technicians to better control contamination during swabbing of critical entry points. These one-step wipes can be used in place of sterile individually wrapped alcohol prep wipes that require tearing of foil/paper pouches in the Primary Engineering Control (PEC) hoods. The low profile peel and reseal cleanroom-compatible packaging provides both quick access and easy closure, and requires minimal space in the PEC. All components of Contec Critical Site Wipes are made in the USA.



Contec Critical Site Wipes' 4" x 4" size make them ideal for use on smaller surfaces where precise application and sterility requirements are mandatory.

Features	Benefits
<b>Presaturated Wipes</b>	<ul style="list-style-type: none"> <li>• Contains consistent amount of 70% USP Grade Isopropanol (IPA)</li> <li>• Allows for process repeatability</li> </ul>
<b>Validated Sterile</b>	<ul style="list-style-type: none"> <li>• Features a Sterility Assurance Level of 10<sup>-6</sup> per AAMI Guidelines</li> <li>• Compliant with USP &lt;797&gt;</li> </ul>
<b>Synthetic, Low Linting Construction</b>	<ul style="list-style-type: none"> <li>• Nonwoven, 100% polyester fabric is cleanroom processed</li> <li>• Cleaner than sterile gauze and other similar products</li> </ul>
<b>Optimal Packaging</b>	<ul style="list-style-type: none"> <li>• Peel and reseal packaging eliminates evaporation concern during storage</li> <li>• Reduces potential risk of contamination from torn foil/paper</li> <li>• 60 wipes per easily recognized pouch is ideal for pharmacy applications</li> </ul>

Part No.	Description	Packaging
HCPS6044	Contec Critical Site Wipes, 4" x 4" (10.2 x 10.2cm) polyester nonwoven wipes, 70% IPA, sterile	60 wipes/pouch, 4 pouches/inner bag, 2 inner bags/outer bag, 3 outer bags/case; 1440 wipes/case